


Fall Detector™

User Instructions (Europe)



D6607016 Issue A Page 1 of 11




Your **fall detector**™ is a battery operated device containing a microprocessor, detector and a radio transmitter. It will automatically trigger an alarm call on your community alarm system if it detects that you have fallen. You can also use it to manually raise an alarm call.

Wearing the fall detector

 Clip the fall detector to your belt or clothing in the upright position shown in the drawing.

It needs to be in the upright position to reduce the possibility of false activation.

Note: An optional Belt and Holster is available (Part Number D6609001C).







Upright position









Incorrect position



Automatic fall detection

- Green light and two "beeps" indicates that the fall detector has been activated.
- The fall detector will then monitor for between 6 and 15 seconds to verify if there could have been a fall. The green light will flash while it does this.
 - If the detector has decided that there is a fall, then a red light and four "beeps" indicates that the *fall detector* is about to trigger an alarm call (because it has detected an impact and has not been upright for between 6 and 15 seconds).
 - If the detector has decided that it was activated by error then the green flashing light will stop after between 6 to 15 Seconds and the unit will revert to normal again.
- The fall detector also features a manual button, which can be used at any time to raise a call for help if required.



Manual operation

- Press the button on the fall detector to raise a manual alarm call.
 - The fall detector will "beep" and a red light indicates the fall detector is triggering an alarm call.
 - If the red light flashes on and off, the fall detector is still triggering an alarm call, but is also indicating to you that the battery needs to be replaced.

D6607016 Issue A Page 4 of 11



- Remember, the fall detector will only detect you falling if you are wearing it.
- The fall detector must be worn at or above the waist, where it will detect the maximum number of falls.



- When removing the fall detector, handle carefully to avoid false calls, place it on an even surface in the upright position and preferably in the stand provided.
- If the fall detector does wake up and beep when you know you have not fallen, you can prevent it raising an alarm by ensuring that it is returned to an upright position.
- It is not recommended that you wear the fall detector whilst in bed, as normal movement in bed will cause the fall detector to trigger.



Battery Replacement

From time to time the battery will need replacing, you will know this because somebody from your monitoring service (if you have one) may inform you, or you may see the red light flashing rather than staying on continuously when a manual alarm call is made.

- The unit automatically tests the battery several times within each 24 hour period.
 It is able to automatically generate an alarm call when low battery is detected.
- If it is necessary for you to change the battery, the fall detector uses a 6 volt Duracell PX28L battery.
- To access the battery, remove the battery cover screw and then the battery cover by sliding it in the direction of the arrow.
- After a year the battery should be replaced. Contact your provider for a replacement. In many cases your provider will change the battery for you. On the back of this user guide is a handy table to



help you remember when the battery should be changed.

- Remove battery and replace it with the new battery, remember to observe the correct orientation. When the battery is first placed in the unit an automatic test call will be generated (no lights will show and no buzzer will sound).
- Replace the sliding cover and the screw.
 Test the fall detector by making a manual call. Inform the operator that you are making a battery test call.
- If the device fails to make a test call then check the orientation of the battery. If this is correct then contact your supplier.



First remove the screw.

Then slide the cover off in the direction of the arrow.

Orientation (copper end).



Do not dispose of the old battery in the bin, please return it to Tunstall (c/o Service Centre) at the address shown overleaf.



Care of the fall detector

- Avoid getting the fall detector wet (if it does then dry it quickly and make a manual test call to confirm it is working correctly).
- If required, the fall detector may be cleaned periodically with a damp cloth and mild detergent.
- It is recommended that you test the fall detector once a month:
 - Make a manual test call.
 - Start with the fall detector vertical.
 Place the fall detector horizontal and then tap within 3 seconds after placing horizontal. After 5 Seconds the fall detector should beep and the green light start to flash. If you put the unit quickly back upright (in its stand, for instance), the alarm call will not be raised.



Notices

This product constitutes a class 2.7 device. It is marked with the CE mark, a notified body number and the 'alert' symbol. These marks have been affixed to demonstrate compliance with the Directive 1999/5/EEC (RTTE Directive).

The 433.15MHz product is intended for use in the following countries:

Germany, Holland, Spain, Denmark, Italy, Sweden, Belgium.

The 469.99MHz product is intended for use only in Germany.

The 173.225MHz product is intended for use in UK, Eire and New Zealand.

The 869.2125MHz product is intended for use in Europe.



D6607016 Issue A Page 9 of 11





01RTTE013B

DECLARATION OF CONFORMITY WITH THE REQUIREMENTS OF EU COUNCIL DIRECTIVE 1999/5/FC.

EQUIPMENT

: Fall Detector 66000/01 (UK 173.225MHz Grey Variant) Fall Detector 66000/02 (UK 173.225MHz Black Variant) Fall Detector 66005/01 (Euro 433.15MHz Grey Variant) Fall Detector 66005/02 (Euro 433.15MHz Black Variant) Fall Detector 66008/01 (Euro 469.99MHz Grey Variant)

Fall Detector 66008/02 (Euro 469.99MHz Black Variant) Fall Detector 67005/49 (Euro 869.2125MHz Grey Variant)

STANDARDS

: EMC : EN50130-4:1995

LVD : EN60950:92 incorporating amendments 1,2, 3, 4 and 11

Radio: EN300-220-3:2000

ASSOCIATED

DOCUMENTS : Summary Information 01RTTE014B

Position

The equipment identified above conforms to the essential requirements of Council Directive 1999/5/EC. The essential radio test suites have also been performed Signed January 2004 Date Technical Virestor

Tunstall Telecom Limited, Whitley Lodge, Whitley Bridge, Yorkshire DN14 0HR

Telephone: 01977 661234 Facsimile: 01977 660562

Internet: www.tunstallgroup.com



Battery Replacement Guide

Battery Inserted Date			Battery Replacement Date		
Day	Month	Year	Day	Month	Year
i					

Remember to change the battery once a year.

Manufactured By

Tunstall Telecom Ltd. Whitley Lodge Whitley Bridge Yorkshire DN14 OHR

Tel: 01977 661234 Fax: 01977 662963 www.tunstallgroup.com **Supplied By**

D6607016 Issue A Page 11 of 11