

DMP version 4.1 Market Release Notes

What's new?

Profile User Level Setting with Time Limitation

It is now possible to set an end date on users access to a specific customer. Settings applied under the Permissions tab on the user.

Edit Isobel Test
✕

OVERVIEW
PERMISSIONS

Permission Profile
Customer Advanced ▼

CUSTOMER NAME	UNTIL	
Tunstall Production	09/12/2020	–
<div> Bristol ▼ </div>	<div> 2020-12-15 </div>	+

DELETE USER
SAVE

You will also notice changes have been made to the appearance of the platform.

If this is the first time you have experienced a DMP upgrade, please note:

- During the upgrade, there will be an expected down time of between 10 and 15 minutes where DMP will be inaccessible.
- Firmware Campaigns and settings changes should not be scheduled to take place within 3 hours of a DMP upgrade.
- Recommended browsers to access DMP are Microsoft Edge, Chrome, and Safari.
- An updated version of the DMP Reference Guide can be found at www.tunstall.co.uk/dmp-guide

For further information please contact the Help Desk on 01977 660204