

# Smart Hub 5.10.4 firmware market release notes

## What's new?

Firmware version 5.10.4 is now available which has the new features listed below.

- Forced cellular network strength scan and selection following every Smart Hub reboot.
- Intruder Alert functionality has been enhanced to provide adequate time to prevent an unintended alarm call when customers return home.

## Smart Hub 4G

We are pleased to announce that **Smart Hub 4G** is also now available, providing compatibility with 4G cellular networks offering the following features:

- VOIP over 4G connection - Voice over data using a non-specific voice channel (for those who cannot access VOLTE)
- VOLTE enabled (network dependent) - Voice over Long-Term Evolution 4G - The 4G network replacing its voice channel with a dedicated 4G voice service.

The device will be preloaded with a new firmware version 5.10.4

5.10.4 firmware is also available to download from DMP for current Smart Hubs (Part No 57200/810.)  
enable them to benefit from the new features.

To

**For further information please contact the Help Desk on 01977 660204**

Tunstall Healthcare Group Ltd, Whitley Lodge, Whitley Bridge, Yorkshire, DN14 0HR  
t: 01977 661234 e: [enquiries@tunstall.com](mailto:enquiries@tunstall.com) w: [tunstall.co.uk](http://tunstall.co.uk) [@TunstallHealth](https://twitter.com/TunstallHealth)

Registered in England No. 06495696